

Thromboresistant¹

ChronoFlex C[®] with BioGUARD[™] Technology

THE **POWERWAND**[™]



**Catheter & Vessel
“Clean and Smooth”**

Control Device



**Mural Thrombus: DVT
“Adherent to the vein for the
length of the device”**

**ChronoFlex C[®] with BioGUARD[™] Technology
is proven in vivo to be thromboresistant¹**

- On the surface of the catheter
- AND**
- On the wall of the blood vessel

Proven to protect where it matters most

Data available on file.

1. POWERWAND is proven in vivo to be thromboresistant with respect to both thrombus on the surface of the catheter and thrombus on the wall of the vein. Based on canine jugular vein thromboresistant study, correlations to clinical applications has not been ascertained.

POWERWAND[™] and BioGUARD[™] are trademarks of Access Scientific, LLC.
ChronoFlex C[®] is a registered trademark of AdvanSource Biomaterials Corp.

access scientific



LC-0586 Rev. A